

# Red Vein Removal

Using Intense Pulsed Light and Laser

a **Fast**  
and effective  
treatment for  
rosacea, facial  
and leg veins

...making light work for YOU!



Before treatment



After one treatment



Before treatment



After one treatment

Courtesy of Dr. Patrick Bowler - Court House Clinic

"After just one treatment the red  
veins on my nose had significantly  
cleared up."

"I have finally lost my thread veins  
and I feel much happier"

## Consultation Procedure

Prior to any treatment the qualified specialist will provide a full consultation. The consultation identifies your expectations and provides a complete explanation of the treatment.

A medical history is taken to confirm suitability for treatment and if you wish to proceed a written informed consent is obtained and a test patch is performed to confirm there are no adverse skin reactions to the light.

**LYNTON**  
LASERS

Tel: +44 (0) 1477 536977  
Fax: +44 (0) 1477 536978  
E-mail: [info@lynton.co.uk](mailto:info@lynton.co.uk)  
Web: [www.lynton.co.uk](http://www.lynton.co.uk)

Ask for more details and book your  
consultation today.

CE  
0120

10<sup>th</sup>  
ANNIVERSARY  
QUALITY • RELIABILITY • INNOVATION

BRINGING YOU  
ADVANCED LIGHT TECHNOLOGY  
IN CLINICAL TREATMENTS

# Red Vein

## Removal

### How does the treatment work?

Intense pulsed light and laser systems release precisely controlled strong pulses of filtered light that are absorbed by the blood vessel, heating it to a point where it is destroyed. Following treatment the vessels quickly clear and are re-absorbed by the body, leaving little or no trace of the original lesion.


### Which types of veins can be treated?

Facial thread veins, Rosacea, Spider Naevi, Port Wine Stains and Leg veins up to 3mm diameter (by laser).

### How is the treatment carried out?

You will be asked to put on goggles to protect your eyes. A layer of chilled gel is applied to the area to guide the light into the skin.

The light applicator is placed on the skin and a short pulse of light is released. The applicator is then moved to the neighbouring area and the process is repeated until the entire area is treated. The chilled gel is removed and the area moisturised.



### Am I suitable for treatment?

The fairer the skin, the easier the treatment. The best response is seen with small thread veins, spider naevi and rosacea. If there are any areas of broken skin, treatment has to be delayed until the area has completely healed.



### How does the treatment feel?

No anaesthesia is required and most patients describe the discomfort as moderate and acceptable.



Sensations vary but most describe it as a very quick hot pin prick.

## ...the light source that works!

### Are there any side effects?

Redness and a warm sensation are normal after treatment. A cool pack may be applied to ease these symptoms. The vessels should fade quickly, and some disappear immediately.

### How long does each treatment take?

Treatments typically take from 15 to 45 minutes, depending on the type and size of lesion.

### How many treatments are necessary?

Typically 1-5 treatments are needed for maximum results, spaced approximately every 4 to 6 weeks.

